



"I have used BloodSTOP iX® in dental extractions where bleeding control was difficult. It costs less than half the price of other products... and it totally dissolves in water so I don't have to worry about a foreign body being left in the surgical area. I have had a 100% success rate with BloodSTOP® iX in over five years of use with hundreds of patients."

Dr. Kenneth Robinson, D.M.D., PC



BloodSTOP iX is the right choice for bleeding control in Dentistry and Oral Surgery

BloodSTOP® iX Absorbent Hemostat is ideally suited to control the significant bleeding associated with tooth extractions and other intra-oral and extra-oral procedures. BloodSTOP® iX is safe, easy to use, and pH neutral. It conforms to the varied spaces and surfaces of the oral cavity, and provides fast-acting bleeding control.

BloodSTOP® iX is a biocompatible, etherized, oxidized, regenerated natural cellulose. This hemostatic matrix rapidly absorbs blood and transforms into a gel to seal the wound with a protective layer, binds and activates platelets, provides support for clot formation, and creates an environment for wound healing.

Uses

- Extractions
- Oral surgery: biopsy, periodontal surgical wounds
- Any other procedures that cause bleeding

Benefits

- Acts in under a minute
 - Is easily inserted into tight oral spaces
 - Odorless, tasteless
 - pH neutral, no added sodium
 - Water soluble for easy removal by irrigation
 - No animal-derived material
 - Biocompatible
 - Safe for all age groups
 - Cost effective
- Sterile, individual use package



"The ability to control bleeding rapidly and effectively in the oral cavity is of major concern to all who perform minor and major intraoral surgical procedures, both oral surgeons and general practitioners. BloodSTOP® iX works very quickly to stop bleeding in a fraction of the time required by traditional cotton gauze and pressure. It is odorless and tasteless, easily applied, and even easier to remove; it washes off with water or saline irrigation."

Dr. George Freedman, DDS



Product Description		Qty per box	Product Code
BloodSTOP Hemostatic Gauze	2"x2"	20/box	BS-10
BloodSTOP iX Absorbent Hemostat	0.5"x2"	24/box	BS-iX-27
BloodSTOP iX Absorbent Hemostat	2"x2"	12/box	BS-iX-14

LifeScience PLUS, Inc.

2520-A Wyandotte Street, Mountain View, CA 94043

Toll Free: 1-877-587-5433 - U.S.A. | Tel: +1-650-565-8172 - International | Fax: 1-650-336-1130

Url: www.lifescienceplus.com. For additional information, contact: sales@lifescienceplus.com

©2012 LifeScience PLUS, Inc., All Rights Reserved. BloodSTOP® is a registered trademark of LifeScience PLUS, Inc.